




# DEFINING A NEW STANDARD IN RADIATION PROTECTION

Say goodbye to heavy lead aprons - the GLiFT™ System is a radiation protection system designed to reduce fatigue, improve ergonomics and provide the protection you deserve.

## Design Features

- GLiFT™ endoskeleton seamlessly integrates with our Standard Vest and Apron
- Revolutionary endoskeleton transfers lead vest weight from your neck, shoulders and back
- Vest and endoskeleton are designed to act as one system — therefore performing like a standard vest with no behavior modification required
- Core material available in lightweight lead and bilayer lead-free alternative
- Multiple colors and print customization available
- Standard 2-piece and 1-piece apron styles available
- Available in multiple lead equivalency weights ranging from 0.25 to .5mm



“

*The weight of the lead is supported by the rigid mantis bars and transferred through the device to the pelvic bone.*

*- Dr. Richard Moore, Jr., MD*

”



# A WEIGHT OFF YOUR SHOULDERS

## (LITERALLY)

“

*I can do a case without it. But, I don't want to.*  
- Interventional Cardiologist, North Carolina

”

### Mantis Arms

Transfer the load of the garment to our proprietary belt system and iliac crest.

### Proprietary Endoskeleton

Engineered for durability, flexibility and individual customization.

### Core Material

Utilizes the highest quality materials available in lead and non-lead options.

### Quality

Crafted from materials worthy of fighter jet design, we hold ourselves to the highest standards possible.

### Comfort

Adjustable belt with memory foam interior and spacer mesh to help regulate user's temperature.



Reduce occupational hazards and long-term negative effects of heavy lead with GLiFT.™

- ▲ Designed by healthcare professionals for healthcare professionals
- ▲ Body contoured design for comfortable all-day wear
- ▲ Economical and proven solution for all healthcare workers



#### TO LEARN MORE:

Email: [sales@theglift.com](mailto:sales@theglift.com) | Phone: 888.765.0570



#### RETURNED GOODS POLICY:

Returned goods authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by contacting GLiFT™ customer service at [sales@theglift.com](mailto:sales@theglift.com)



GLiFT™ is a registered trademark of Salus Scientific, Corp.



**SALUS  
SCIENTIFIC**