

The patented GLiFT™ is designed to provide best-in-class radiation protection and help "carry the load" for healthcare professionals.



The GLiFT™ Difference

- Revolutionary design transfers lead weight to hips
- Designed by healthcare professionals
- ✓ Body contoured design for comfortable all-day wear
- Economical and proven solution for all healthcare workers





Design Features

- GLiFT™ endoskeleton seamlessly integrates with our Standard Vest and Apron
- Revolutionary endoskeleton transfers lead vest weight from your neck, shoulders and back
- Vest and endoskeleton are designed to act as one system — therefore performing like a standard vest with no behavior modification required
- Designed with extensive input from practicing healthcare professionals
- Core material available in lightweight lead and bilayer lead-free alternative
- Multiple colors and print customization available
- Standard 2-piece and 1-piece apron styles available
- Available in multiple lead equivalency weights ranging from 0.25 to 1.0 mm



Your Back Will Thank You!

www.theGLiFT.com



GLIFT Part Number **Description** Radiation Protection Vest that **VST-SM** Incorporates a gravity longitudinal **Vest Part VST-MD** individualized force transfer system. Number **VST-LG** Available in lead and lead-free core **VST-XL** material. • Radiation protection kilt. Available in **Standard** KLT-SM lead and lead-free core material. Vest / Kilt **Kilt Part KLT-MD** System KLT-LG Number KLT-XL Radiation protection thyroid shield. TS-SM Available in lead and lead-free core **Thyroid** TS-MD Shield material. TS-LG TS-XL Radiation protection system that OPA-SM incorporates a gravity longitudinal Apron OPA-MD individualized force transfer system. Part OPA-LG Available in lead and lead-free core Number OPA-XL **One-Piece** material. Apron **System** • Radiation protection thyroid shield. TS-SM Available in lead and lead-free core **Thyroid** TS-MD Shield TS-LG material.

Protected by U.S. Patent No. 10,729,195 and other U.S. and foreign patents pending.

TS-XL

Latex and DEHP Status: Latex, BPA, and DEPH Free

REMIT TO:

Salus Scientific Corporation 835 Coleman Blvd., Suite 101 Mount Pleasant, SC 29464

CONTACT US:

Email: sales@theglift.com Phone: 888.765.0570



88-1543518

F RETURNED GOODS POLICY:

Returned goods, authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by contacting GLiFT™ customer service at sales@theglift.com





